

Focus Biomolecules Protease Inhibitors

Target	Product Name	Catalog #	CAS	Size A	Size B
ACE	Captopril	10-2253	62571-86-2	1 g	5 g
Aminopeptidase	Puromycin 2HCl	10-2100	58-58-2	50 mg	250 mg
Aminopeptidases	Bestatin HCl	10-1286	202475-60-3	5mg	25 mg
Aspartic proteases	Pepstatin A	10-1368	26305-03-3	5 mg	25 mg
BACE1	Asiatic acid	10-3020	464-92-6	10 mg	50 mg
Calpain	Ac-Leu-Leu-Methional (ALLM)	10-1349	136632-32-1	5 mg	25 mg
Calpain	PD-150606	10-1515	179528-45-1	5 mg	25 mg
Calpain	Ac-LLnL-CHO	10-3024	110044-82-1	5 mg	25 mg
Calpain	MDL-28170	10-3005	88191-84-8	5 mg	25 mg
Caspase (pan)	ZVAD-FMK (OMe)	10-2614	187389-52-2	1 mg	
Chymotrypsin	SL-01 (Z-Phe-CMK)	10-1326	26049-94-5	50 mg	250 mg
Chymotrypsin-like Ser proteases	TPCK	10-2290	402-71-1	100 mg	500 mg
Cruzin	Bromo-PTC	10-1433	NF	10 mg	50 mg
Cysteine proteases	E-64	10-1345	66701-25-5	5 mg	25 mg
Cysteine proteases	E-64d	10-1347	88321-09-9	1 mg	5 mg
DPP-4	P32/98	10-1142	136259-20-6	10 mg	50 mg
DPP-4	Vildagliptin	10-2887	274901-16-5	10 mg	50 mg
DUB	b-AP15 (NSC-687852)	10-1340	1009817-63-3	5 mg	25 mg
DUB	LDN-57444	10-1303	668467-91-2	10 mg	50 mg
DUB	NSC-632839	10-1385	157654-67-6	10 mg	50 mg
DUB	PR-619	10-1485	2645-32-1	5 mg	25 mg
DUB	TCID	10-2428	30675-13-9	10 mg	50 mg
DUB	WP-1130 (Degrasyn)	10-1351	856243-80-6	5 mg	25 mg
DUB (USP-14)	IU1	10-1197	314245-33-5	10 mg	50 mg
DUB (USP-7)	P5091	10-1422	882257-11-6	5 mg	25 mg
Elastase	Sivelestat	10-1102	127373-66-4	5 mg	25 mg
Elastase	Elastatinal	10-2818	51798-45-9	5 mg	25 mg
g-Secretase	DAPT	10-4514	208255-80-5	5 mg	25 mg
g-Secretase	JLK-6	10-1568	62252-26-0	5 mg	25 mg
MetAP-2	Fumagillin	10-1146	23110-15-8	1 mg	5 mg
MMP	Ilomastat (GM-6001)	10-1339	142880-36-2	1 mg	5 mg
MMP-13	CL-82198	10-1028	307002-71-7	5 mg	25 mg
MMP-14	NSC-405020	10-1462	7497-07-6	5 mg	25 mg
MMPs	NNGH	10-1045	161314-17-6	5mg	25 mg
Neutral endopeptidase	Phosphoramidon 2Na	10-2548	119942-99-3	5 mg	25 mg
Proteasome	Bortezomib	10-2120	179324-69-7	5 mg	25 mg
Proteasome	Carfilzomib	10-1477	868540-17-4	1 mg	5 mg
Proteasome	Lactacystin	10-1314	133343-34-7	200 ug	1 mg
Proteasome	MG-132	10-1309	133407-82-6	5 mg	25 mg
Ser, Cys & Thr proteases	Leupeptin hemisulfate	10-1346	103476-89-7	25 mg	100 mg
Serine and cysteine proteases	Antipain 2HCl	10-2734	37682-72-7	5 mg	25 mg
Serine proteases	AEBSF HCl	10-1348	30827-99-7	50 mg	
Serine proteases	Chymostatin	10-2722	9076-44-2	5 mg	25 mg
Serine proteases	Nafamostat Mesylate	10-1089	82956-11-4	10 mg	50 mg
Thrombin	Argatroban	10-2809	74863-84-6	10 mg	50 mg
Trypsin & serine proteases	Benzamidine HCl	10-2687	1670-14-0	5 g	25 g
Zn metalloproteases	1,10-Phenanthroline HCl	10-3017	3829-86-5	1 g	5 g

Available in Europe at...

tebu-bio.com
Innovative Lab Service & Reagents